Dräger DrugTest® 5000 Test-Kit

REF 83 19 830

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en Instructions for Use

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Directive 98/79/EC

Dräger Safety AG & Co. KGaA
Revalstraße 1
D-23560 Lübeck
Germany
Tel. +49 451 8 82-0
Fax +49 451 8 82-20 80
www.draeger.com

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en – Instructions for Use

For Your Safety
Observe the instructions for use
Any use of the Dräger DrugTest® 5000 Test Kit requires precise knowledge and observation of the instructions for use.
The components of the test kit are only intended for the usage described.
Dräger Safety AG & Co. KGaA

Intended Use

The Dräger DrugTest® 5000 Test System consists of the Dräger DrugTest® 5000 Analyzer and the Dräger DrugTest® 5000 Test Kit. The test system is designed for the simultaneous, qualitative detection of substances and substance classes in human saliva for diagnostic purposes (in-vitro diagnostic) or for forensic application.
The Test Kit 6 Panel (Ref. 83 19 830) is designed for the qualitative detection of cocaine, opiates, benzodiazepines, THC, amphetamine and methamphetamine, in cooperation with the DrugTest 5000 Analyzer.
The Dräger DrugTest 5000 Test System is a qualitative measurement procedure. It detects the presence of specified substances or their metabolites in the sample, below or under a cut-off, and therefore only provides a preliminary analytical result (screening procedure). To obtain a confirmed analytical result a more specific method must be used. As a rule the preferred method is gas chromatography/mass spectrometry (GC/MS). All results from the DrugTest 5000 Test System require a professional analysis above a cut-off concentration. Therefore, it only provides a preliminary analytical result (screening procedure). To obtain a confirmed analytical result a more specific method must be used. As a rule the preferred method is gas chromatography/mass spectrometry (GC/MS).
All results from the DrugTest 5000 Test System require a professional analysis above a cut-off concentration. Therefore, it only provides a preliminary analytical result (screening procedure). To obtain a confirmed analytical result a more specific method must be used. As a rule the preferred method is gas chromatography/mass spectrometry (GC/MS).
For details of how to use the analyzer, see the Instructions for Use for the DrugTest 5000 Analyzer.

Dräger DrugTest® 5000 Cut-Off

Drug	Calibrator	ng/ml
COC	Cocaine	20
OPI	Opiate	20
BENZO	Benzodiazepine	15
THC	Delta-9-Tetrahydrocannabinol	5
AMP	Amphetamine	50
MAMP	Methamphetamine	35

Warnings and precautionary measures

The packaging with the test cassette must be completely sealed. Do not use the test cassette if the seals are not intact.
To prevent the contamination of the oral fluid collector, only open the packaging immediately before use.
Do not touch the oral fluid collector with bare fingers before or after specimen collection and observe standard hygiene procedures. The Test Kit can only be processed using the Dräger DrugTest® 5000 Analyzer!

Test principle

The Dräger DrugTest® 5000 Test System is based on the immuno assay principle of competitive displacement. Drugs found in the saliva compete with drugs on the test membrane to bind onto microparticles transferred by antibodies.

Samples: The Test Kit is designed to be used for samples collected with the oral fluid collector. The sample must not be treated in a special manner. The saliva is absorbed directly into the porous collector, which is integrated in the test cassette. Once the samples have been taken, the analysis is started by inserting both the test cassette and the cartridge into the analyzer. The analyzer automatically transfers the sample to the test cassette, whereby the test development is triggered.

Test: The saliva sample interacts with microparticles (transferred with antibodies) and with drug conjugates on the test membrane. If the sample is free of drugs, the antibodies can react freely with the drug conjugates, resulting in a signal being generated on the test membrane. If drug exist in the sample, they bind onto the microparticles transferred with the antibodies, whereby the signal generated becomes weaker. The analyzer detects the signal generated by a specific sample, and decides whether this signal comes from a positive preliminary sample ("not negative").

Quality control: Another antibody/antigen reaction, which is independent of the sample, is integrated in each test membrane. If the sample processing is successful, antibodies bind onto the test membrane and the antigens bind onto the microparticles, thus generating a test signal. This signal is also detected by the analyzer and is used to decide whether or not a test is valid.

Data management: After analysis, test results are shown in the analyzer's display and are stored in the analyzer's memory. All data can be printed. For further details on the analyzer's data management options, see the technical manual for the Dräger DrugTest 5000 Test System¹⁾.

Storage and stability

The cartridges must be stored in their original pouches at between 4 and 30 °C (39 and 86 °F). Do not use the cartridges if the pouches are damaged (e.g. if there is a hole or a tear).
The cartridge must be used immediately after opening the pouches.
Do not use the test kit once the expiry date has passed.

Scope of delivery

20 individually packed test cassettes
20 cartridges
1 Instructions for Use

Material required but not included

Dräger DrugTest® 5000 Analyzer (order No. 83 19 900) for processing and analysing the test kit.
Additional, recommended material
Dräger DCD 5000 (order No. 83 19 910),
Disposable gloves (e.g. latex or nitrile gloves).
For technical support, contact your local IVD dealer, Dräger Safety or Technical Support.

Quality control

Quality control testing at regular intervals is good analytical practice and may be required by federal, state, or local guidelines. Always check with the appropriate licensing or accrediting bodies to ensure your quality control program meets the established standards.

Internal control: An integrated process control is performed for each analysis (see also section "Test principle"). If a test was developed correctly, a control signal is generated. If this signal is not generated, the analyzer does not evaluate the corresponding analyses, but issues the result of this analysis as "invalid".
External control: The characteristics of the test mean that a quality control can be performed on site. A drug-free sample is taken and processed. The result for all analyses with this sample must be NEGATIVE. For information on positive commercial controls for this test, please contact the Dräger Safety technical support.

Test execution

Preparation

- Inform the test subject not to consume anything ten minutes prior to taking the sample (e.g. food, drink, tobacco, chewing gum).
- Ensure the specified ambient conditions: ambient air temperature between 5°C and 40°C; relative humidity between 5% and 95% r.h.
- Ensure that the test cassettes in the pouches and cartridges are at ambient temperature (if necessary, wait 15 min for temperature compensation).
- Only open pouches immediately before use, remove test cassette from the pouch.

Collecting the oral fluid sample

- Remove the protective cap from the oral fluid collector of the test cassette and hand the test cassette to the subject – figure 1.
- Instruct the subject to place the oral fluid collector inside the cheeks and move it carefully from one side of the mouth to the other – figure 2. Do not chew or suck on the oral fluid collector.
Observe the collection!
– The sample quantity display turns blue when a sufficiently large amount of sample has been taken. The subject can then remove the oral fluid collector from the mouth – figure 3.
– If the sample is not analyzed immediately after collection, you must read the notes in the chapter "Batch mode".

Performing the evaluation

- Take the test cassette from the subject, open the door of the Analyzer and insert the test cassette into the lower compartment, until it engages audibly – figure 4.
- Insert the cartridge into the upper compartment of the Analyzer, until it engages audibly – figure 5.
- Close the door of the analyzer. The analyzer then starts the analysis automatically.
– Once the analysis is complete, the analyzer displays the results of each tested substance on the screen. For more details, please see the Instructions for Use for the analyzer.
- After successful analysis, remove the test cassette with the cartridge attached from the analyzer and dispose of the used test kit – figure 6.

Batch mode

If the analysis of the oral fluid sample is not to be performed directly after taking the sample, the test cassette can be stored at room temperature for up to 4 hours (evidence of THC) or 8 hours (evidence of all other drugs). In this case the protective cap must be placed back on the oral fluid collector. The protective cap must be removed from the oral fluid collector again before analysis. Execute the analysis as described under "Executing analysis".
Caution: Failure to adhere to the maximum holding time can result in invalid or incorrect results.

Confirmation of results

The Dräger DrugTest® 5000 Test System is a qualitative measurement procedure for detecting specific substances or their metabolites in the sample above a cut-off concentration. Therefore, it only provides a preliminary analytical result (screening procedure). To obtain a confirmed analytical result a more specific method must be used. As a rule the preferred method is gas chromatography/mass spectrometry (GC/MS).
To obtain a confirmed result, it is necessary for a second sample to be collected and analysed by a recognised laboratory. The second sample is collected using the Dräger DCD 5000 (order no. 83 19 910).

Limitations

- The Dräger DrugTest® 5000 Test Kit has been designed for analysing oral fluid samples taken using the test kit's oral fluid collector. The use of other liquids (e.g. urine, sweat, buffered calibrators) or oral fluid samples collected in a different manner can lead to inaccurate results. Only use accessories permitted by Dräger for the Dräger DrugTest® 5000 Test System.
- A positive test result only shows that the concentration of the substance lies over the respective test kit limit value for the substance. A negative result can show that the sample is either free from drugs or that the concentration lies below the limit value.
- The reagents used in this test are specifically designed to detect the drugs specified. Nevertheless, structurally similar prescribed and over-the-counter drugs can react with the antibody reagents and cause false positive results. The characteristics of immuno assay mean that cross-reactive substances or unspesific interactions can disrupt the test and lead to inaccurate results. Technical and procedural errors can also contribute to incorrect results.
- In rare cases, inaccurate results can also result from unspesific interactions (physiological deviation, state of health, contamination of the saliva through food or other substances) (see also "Disruptive substances"). Technical and procedural errors can also contribute to incorrect results.
- The extraction of an oral fluid sample can be monitored precisely. As a result, specimen adulteration is unlikely. However, if you suspect that the sample may have been adulterated, repeat the test with a new test kit.

Disruptive substances

According to our tests, the following formulations do not lead to malfunctions of the Dräger DrugTest® 5000 analyser if ten minutes have passed since they were consumed: Cola, toothpaste, American cranberry juice, aseptic mouthwash, orange juice, water, and cough syrup.
The substances listed in table 1 (T1) were tested in quantities of 10,000ng/ml, and do not influence the Dräger DrugTest® 5000 analyses.

Specificity

Substances similar to drugs that can be detected by the system were investigated to determine possible cross-reactivity. The quantities listed in table 2 (T2) represent maximum concentration limits for the named substances in order to generate a positive test result.

Accuracy

In total, 100 oral fluid samples were analysed using the Dräger DrugTest® 5000 Test System. A second sample was taken at the same time and analysed by GC/MS. The clinical performance of the Dräger DrugTest® 5000 Test System is summarised in table 3 (T3).

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T2 – Specificity	
Cocaine related compounds	[ng/ml]
Benzoylcegonine	200
Cocacethylen	500
Ecgonine methyl ester	>10000

Benzodiazepine related compounds	
7-Aminoflurazepam	[ng/ml]
Alprazolam	250
Bromazepam	150
Chlordiazepoxide	25000
Clozapepam	90
Desalkylflurazepam	45
Flunitrazepam	20
Flurazepam	30000
Lorazepam	200
Midazolam	100
Nitrazepam	30
Nordiazepam	45
Oxazepam	50
Przepam	100000
Temazepam	15

THC related compounds	
11-Nor-9-Carboxy-9-THC	[ng/ml]
Cannabinol	100000
Cannabinol	100000

Amphetamine related compounds	
MDA	[ng/ml]
MDA	100
MDMA	>10000
MDMA	>10000
MBDB	>10000
Phentermine	500
Si(+)-Methamphetamine	30000
Tyramine	5000

Methamphetamine related compounds	
MDA	[ng/ml]
MDA	1000
MDMA	>10000
MDMA	75
MBDB	35
Si(+)-Amphetamine	30000
Procain	5000
Pseudoephedrine	10000

T3 – Clinical Performance			
Cocaine			
Cocaine-Equivalents			
Dräger DrugTest® 5000	GC/MS [ng/ml]		
	<8 ng	8 – 20 ng	>20 ng
positive	0	1	20
negative	95	5	0
Overall agreement between Dräger DrugTest® 5000 and GC/MS: 98 % (Cut-off GC-MS: 20 ng/ml)			

Opiates			
Morphin-Equivalents			
Dräger DrugTest® 5000	GC/MS [ng/ml]		
	<8 ng	8 – 20 ng	>20 ng
positive	0	7	27
negative	86	5	1
Overall agreement between Dräger DrugTest® 5000 and GC/MS: 94 % (Cut-off LC-MS/MS: 10 ng/ml)			

Benzodiazepine			
Diazepam-Equivalents			
Dräger DrugTest® 5000	HPLC/MS [ng/ml]		
	<4 ng	4 – 15 ng	>15 ng
positive	0	3	17
negative	84	3	2
Overall agreement between Dräger DrugTest® 5000 and LC-MS/MS: 97 % (Cut-off LC-MS/MS: 7.5 ng/ml)			

Tetrahydrocannabinol			
Delta-9-Tetrahydrocannabinol			
Dräger DrugTest® 5000	GC/MS [ng/ml]		
	<5 ng	5 – 15 ng	>15 ng
positive	0	2	23
negative	84	4	2
Overall agreement between Dräger DrugTest® 5000 and GC/MS: 94 % (Cut-off GC-MS: 2.5 ng/ml)			

Amphetamine			
Amphetamine-Equivalents			
Dräger DrugTest® 5000	GC/MS [ng/ml] *		
	<20 ng	20 – 50 ng	>50 ng
positive	1	1	19
negative	94	4	0
Overall agreement between Dräger DrugTest® 5000 and GC/MS: 98 % (Cut-off GC-MS: 50 ng/ml)			

Methamphetamine			
Methamphetamine-Equivalents			
Dräger DrugTest® 5000	GC/MS [ng/ml] *		
	<14 ng		

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© Dräger Safety AG & Co. KGaA